

Application Note

Sandwich ELISA for Detection of Collagenase

Protocol for sandwich ELISA suitable for quantification of Collagenase NB 4 Standard Grade, Collagenase NB 5 Sterile Grade or Collagenase NB 6 GMP Grade in solutions. Concentrations as low as 100 ng/ml can be detected with antibodies used in this protocol.

Equipment:

- > pH meter
- > Multichannel pipette with tips
- > Microtitre plate (transparent, 96 wells) with lid
- > Microplate reader

Reagents:

- > Collagenase standard: Collagenase 4 Standard Grade, Collagenase NB 5 Sterile Grade, or Collagenase NB 6 GMP Grade
- > Capture antibody: Anti-Collagenase (*Clostridium histolyticum*), polyclonal antibody from sheep
- > Tracer antibody: Anti-Collagenase (*Clostridium histolyticum*), polyclonal antibody from sheep, conjugated with horseradish peroxidase

Buffers and solutions:

- > Coating buffer: 50 mM Na₂CO₃, pH 9.6 (adjust with 1 M HCl)
- > Phosphate buffered saline (PBS), 1 x, pH 7.4 (adjust with 0.1 M NaOH)
- > Wash buffer: PBS with 0.05 % (v/v) Tween 20
- > Blocking solution: PBS with 0.1 % bovine serum albumin
- > Substrate solution:
0.2 M Na₂HPO₄ and 0.1 M citric acid in H₂O, pH 5.0 (adjust with 0.1 M NaOH)
with 0.04 % (w/v) ortho-phenyldiamine dihydrochloride
Attention: ortho-phenyldiamine dihydrochloride is light-sensitive!
Aliquots of the substrate solution can be stored at -20 °C for six months.
- > Hydrogen peroxide solution 30 % (w/w) in H₂O
- > Stop solution: 2 M H₂SO₄ in H₂O

For information and samples please contact the Nordmark Biochemicals Team

Phone: +49 4122 712 560 | Fax: +49 4122 712 286
E-Mail: collagenase@nordmark-biochemicals.com
Pinnauallee 4 | 25436 Uetersen | Germany
www.nordmark-biochemicals.com

Ordering

Nordmark Biochemicals
Phone: +49 4122 712 560
Fax: +49 4122 712 286
E-Mail: order@nordmark-biochemicals.com



Application Note

Sandwich ELISA for Detection of Collagenase

Procedure:

1. Dilute the capture antibody in coating buffer to a final solution of 10 µg/ml. Coat the wells by adding 50 µl (500 ng capture antibody) into each well.
2. Seal the microtitre plate with a lid and incubate for at least 16 h at +2 to +8 °C.
3. Invert the plate to empty and gently tap out the residual liquid on a paper towel. Add 300 µl wash buffer per well and incubate for 30 s at room temperature. Wash three times in total. After the last washing step gently tap out residual liquid on a paper towel.
4. Add 300 µl blocking solution per well.
5. Seal the plate and incubate for 30 min at 37 °C.
6. Wash the plate three times (see step 3).
7. For generating a standard curve prepare on ice dilutions of collagenase standard (in triplicate) from ~1 ng/ml to 10 µg/ml in fresh buffer which was prepared for washing of the cells. If this buffer is not available, use wash buffer [PBS with 0.05 % (v/v) Tween 20] for dilution of standards. Include buffer without collagenase as blank. Dilute on ice samples in the same buffer which was used for dilution of standards.
8. Add 100 µl of above mentioned solutions per well.
9. Seal the plate and incubate for 60 min at 37 °C and 300 rpm (orbital shaker).
10. Wash the plate three times (see step 3).
11. Dilute the tracer antibody in washing buffer to a final solution of 5 µg/ml. Add 50 µl (250 ng tracer antibody) per well.
12. Seal the plate and incubate for 60 min at 37°C and 300 rpm (orbital shaker).
13. Wash the plate three times (see step 3).
14. Add 1 µl H₂O₂ per 1 ml substrate solution, mix and add 100 µl per well.
15. Seal the plate and incubate for 20-30 min at room temperature in the dark.
16. Stop reaction by adding 50 µl stop solution per well.
17. Measure the optical density at 492 nm with a microplate reader within 30 min after adding stop solution. Plot standard curve and use it to quantify the collagenase in the buffer used for washing of the cells.

Application Note

Sandwich ELISA for Detection of Collagenase

Ordering Information

Product	Cat. No.	Pack size
Anti-Collagenase (C. histolyticum), polyclonal antibody from sheep	S5805001	100 µg
Anti-Collagenase (C. histolyticum), polyclonal antibody from sheep, conjugated with horseradish peroxidase	S5805101	100 µg
Collagenase NB 4 Standard Grade	S1745402	500 mg
Collagenase NB 4 Standard Grade	S1745401	1 g
Collagenase NB 4 Standard Grade	S1745403	5 g
Collagenase NB 5 Sterile Grade	N0002778	1 g
Collagenase NB 6 GMP Grade	N0002880	100 mg
Collagenase NB 6 GMP Grade	N0002779	1 g

The user of this protocol is solely responsible and liable.