New Product

Need some extra Power that's extra Gentle?

Try Collagenase NB 7D Gentle Grade for



3D model of lymphoma cell

Sensitive Cells


Dendritic cells, diverse tumor cells and others

• Preservation of Receptor Functionality Low tryptic activities for sensitive conditions

Research Applications

Affordable pricing with excellent quality

The naturally-expressed enzyme mixture of collagenases can be modified within the manufacturing process to create products adapted for different applications. Now, Nordmark provides a complete spectrum of collagenase products to match your needs.



Collagenase NB 7D Gentle Grade is a partially purified collagenase with slightly **higher collagenase activity but low tryptic activities.** This allows preservation of cell-surface proteins or other delicate structures. As a result, Collagenase NB 7D Gentle Grade is suitable for isolation of a broad variety of cells from human or animal tissues where extra care is required.



Product Information

| Enzyme | Collagenase NB 7D Gentle Grade | | |
|-------------------|--|------------------------------|--|
| Origin | Clostridium histolyticum | | |
| Composition | Contains class I and class II collagenases and reduced trypsin-like activities | | |
| Molecular weight | 70 – 120 kDa (collagenases) | | |
| Special features | Slightly higher collagenase activity with reduced trypsin-like activities | | |
| Character | White to beige powder, lyophilisate | | |
| Activities | Collagenase (PZ activity acc. to Wünsch, 25 °C) | \geq 0.2 U/mg lyophilisate | |
| | Trypsin-like activities (BAEE) | ≤ 0.2 U/mg lyophilisate | |
| | Ratio Collagenase : Trypsin | ≥ 3.0 | |
| Long term storage | +2 to +8 °C | | |

Collagenase NB 7D Gentle Grade can be used for a variety of applications where gentle conditions are beneficial:

- Tumor cells
- Lung cells
- Lymphocytes
- Endothelial cells Dendritic cells
 Splenocytes
- Leukocytes
- Many more...

Contact us for more information or for free samples at:

collagenase@nordmark-biochemicals.com



Ordering Information

| Enzyme | Cat. No. | Pack Size |
|--------------------------------|----------|-----------|
| Collagenase NB 7D Gentle Grade | S1747401 | 500 mg |
| | S1747402 | 1 g |
| | S1747403 | 5 g |

For information and samples please contact the Nordmark Biochemicals Team

Phone: +49 4122 712 560 Fax: +49 4122 712 286 E-Mail: collagenase@nordmark-biochemicals.com Pinnauallee 4 25436 Uetersen Germany www.nordmark-biochemicals.com

Ordering

Nordmark Biochemicals Phone: +49 4122 712 560 Fax: +49 4122 712 286 E-Mail: order@nordmark-biochemicals.com



Nordmark Pharma GmbH | Pinnauallee 4 | 25436 Uetersen | Germany | Phone: +49 4122 712 0 | info@nordmark-pharma.de